

GrainPro® Hermetic Bag Premium™

(ZIPPER M-L)



1

Fold the GrainPro Hermetic Bag lengthwise and insert it inside the outer bag.

NOTE: Refer to suggested outer bag dimensions below.

2

Unfold the GrainPro Hermetic Bag and fold the extra liner to prevent damage.

3

Pour the dry commodity until it reaches the MAXIMUM FILL LINE.

Note: Overfilling can cause the zipper to reopen.

4

Unroll and remove air by pressing down on the commodity.

OPTIONAL: Add an oxygen absorber to lower the oxygen level immediately.

If storing pulses, 1 sachet of oxygen absorber per 50 kg of the commodity must be added.

5

Press one end of the track and push the slider from the top end of the zipper. To seal, apply slight pressure and move the slider to the other end.

6

Fold down the zipper twice for better sealing.

7

Close the outer bag by tying or sewing at a 3 cm distance from the "zipped" inner liner.

8

Stick the 'No Hooks' label on the outer bag.

9

For future use, remember to clean, dry, and fold the used and undamaged GrainPro Hermetic Bag.

10

Sufficiently low permeation against gases, moisture and water.

11

No need to use fumigants.

**For a Foodful Future
GrainPro: The Best Choice for Quality Preservation**

For more info, scan the QR code or email us at customercare@grainpro.com

Outer Bag recommended dimensions:
 GrainPro® Hermetic Bag Premium™ Zipper Medium, cm (inch) 65-70 x 64-65 (25.6-27.5 x 25.2-25.6)
 GrainPro® Hermetic Bag Premium™ Zipper Large, cm (inch) 80-85 x 70-72 (31.5-33.4 x 27.5-28.3)

GRAINPRO® HERMETIC BAG PREMIUM™ ZIPPER		
SIZE	Approximate Capacity kg (lbs) Based on wheat	Dimensions, L x W, cm (i nch)
Medium	30 (66)	80 x 65 (32 x 26)
Large	69 (152)	110 x 72 (43 x 28)

Disclaimer: The approximate capacity provided above are based upon our computation. Variance upon actual usage may apply.