GrainPro® Hermetic Bag Premium™



(ZIPPER XS-S)



Unfold the GrainPro Hermetic Bag Premium and fold the extra portion to prevent damage.



Pour in the dry grains. Leave 10 cm of the liner free for proper sealing.



Unroll and remove air by pressing down on the grains

OPTIONAL: Add an oxygen absorber to lower the oxygen level immediately.

If storing pulses, 1 sachet of oxygen absorber per 50 kg of the commodity must be added.



Press one end of the track and push the slider from the top end of the zipper. To seal, apply slight pressure and move the slider to the other end.

WARNING: Do not leave the zipper slider attached to the bag. It must be removed and kept for future use.




Fold the zipper down twice for better sealing.



Insert the folded portion onto the side of the outer box.



For future use, remember to clean, dry, and fold the used and undamaged GrainPro Hermetic Bag.



No need to use fumigants.

Sufficiently low permeation against gases, moisture and water.



For a Foodful Future GrainPro: The Best Choice for Quality Preservation

For more info, scan the QR code or email us at customercare@grainpro.com

70 x 40 (28 x 16)

Outer Bag recommended dimensions:

Hermetic[®] Bag Premium[™] Zipper Extra Small, Designed as a liner for cardboard box Hermetic[®] Bag Premium[™] Zipper Small, Designed as a liner for cardboard box

GRAINPRO® HERMETIC BAG PREMIUM™ ZIPPER Approximate Capacity kg (lbs) Based on green coffee 5 (11) Dimensions, L x W, cm (inch) 50 x 30 (20 x 12)

Disclaimer: The approximate capacity provided above are based upon our computation. Variance upon actual usage may apply.

15 (33)

Small

SIZE

Extra Small